






Note: Always swab self sealing port with sterile alcohol prior to accessing with a sterile syringe

<p>Step 1:</p>  <p>Draw 6mL of Sodium Citrate Anticoagulant into 60mL Syringe</p>	<p>Step 2:</p>  <p>Draw 54mL whole blood from the patient, filling the syringe to 60mL</p>	<p>Step 3:</p>  <p>Load anticoagulated whole blood into the Concentrating Device</p>	<p>Step 4:</p> <p>Counterbalance and process the Concentrating Device at</p> <p>Executive Series Centrifuge 5 minutes 4400 RPM</p> <p>Platinum Series Centrifuge Pre-Programmed Setting: "AbsolutePRP/BMC"</p>
<p>Step 5:</p>  <p>Piston touches the RBC Interface (trapping plasma & platelet buffycoat inside)</p> <p>Attach the 60mL syringe and aspirate the plasma until the piston touches the RBC interface (trapping the plasma & platelet buffycoat inside), then stop aspirating</p>		<p>Step 6:</p>  <p>7mL PRP (Plasma & Platelet buffycoat)</p> <p>Attach the 12mL syringe and aspirate 7mL PRP (Rotate the syringe to re-suspend the platelet buffycoat into the plasma)</p>	